

Facility/Laboratory Name: _____

INCUBATOR TEMPERATURE LOG

OFFICIAL RECORDS FOR THE PERIOD _____ TO _____

Make/ Model: _____ Serial Number: _____ Unit Id#/ Identification#: _____ Temperature Range Used : _____

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM Top																															
AM Bottom																															
PM Top																															
PM Bottom																															

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM Top																															
AM Bottom																															
PM Top																															
PM Bottom																															

Month	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM Top																															
AM Bottom																															
PM Top																															
PM Bottom																															

1. Size adequate for work load
 2. Maintains appropriate temperature range for plate incubation.
 3. Not to be used to store food or drink for consumption.
 4. Record temperature (corrected) daily top and bottom shelves if applicable, in AM and PM, thermometer with bulb immersed in antifreeze liquid (in sealed containers)
 5. Record temperature daily to the nearest 0.1°C, maintain records. Note any corrective actions taken.
- Example:

32.5°C MH
