

Facility/Laboratory Name: \_\_\_\_\_ **FREEZER TEMPERATURE LOG**

**OFFICIAL RECORDS FOR THE PERIOD** \_\_\_\_\_ **TO** \_\_\_\_\_

Make/ Model: \_\_\_\_\_ Serial Number: \_\_\_\_\_ Unit Id#/ Identification#: \_\_\_\_\_ Temperature Range Used : \_\_\_\_\_

**Month**

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM																															
PM																															

**Month**

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM																															
PM																															

**Month**

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
AM																															
PM																															

1. Size adequate for work load
2. Maintains -15°C or below
3. Used for storage of frozen milk products, controls, media and reagents only. Not to be used to store food or drink for consumption.
4. Record temperature (corrected) daily, in AM and PM, thermometer with bulb immersed in antifreeze liquid ( in sealed containers)
5. Record temperature daily to the nearest 0.1°C, maintain records. Note any corrective actions taken.

Example:

-17.5°C
MH